

PRO₂XY[®] VAC

VPSA Technology : Vacuum Pressure Swing Adsorption



Oil-less low pressure compressor to feed the oxygen generator.
No need of air treatment system as for PSA.
No oil pollution.

Set of pneumatic valves developed internally and extremely reliable.



Low maintenance cost due to 100% oil free process from 0,5 bar to 12 bar.

Better Total Cost Ownership comparing to PSA

Pressure/vacuum tanks made in our internal boiler and filled with zeolite lithium basis.



VPSA oxygen generators

Oil-less vacuum pump to regenerate the zeolite and participate to energy saving.



Built-in electrical connections on each machine for plug and play installation.



Oil-less scroll oxygen compressor from 0,5 bar to 6 bar driven by variable frequency converter for energy saving.



Oil-less piston booster from 6 to 12 bar for double stage pressure pipeline.



Our new VPSA generators will give you a technological edge for oxygen on-site production.

VPSA is the acronym for Vacuum Pressure Swing Adsorption.

Vacuum : To help the regeneration / desorption of the molecular sieve
Pressure : To feed the generator molecular sieve
Swing : The steady changeover of production bed
Adsorption : Is based on the ability for porous materials like zeolites to bind gases through their large surface areas.

But it also stands for Innovation :

- No more energy waste to compress nitrogen at several bars.
- Special zeolite with high affinity for water vapor, removing the need of dryer.
- No oil : vacuum pump and blower are oil-less rotary lobe units.

Mil's has drawn from its experience of PSA generators to develop this new concept for reliable production on the hospital site, independent from cylinder or liquid supply.

The adsorption-desorption phenomena described on the next page is similar to the PSA one with its 2 beds systems. However the cycle being carried out at lower pressures, the power consumption will decrease of more than 60% for a constant production of oxygen at 95%.

The new Procom 3 touchscreen panel ensures a trouble free production, with control and record of all the parameters. Paired with the Vigifluid remote monitoring, every connected user can follow and trace back the oxygen quality and events. The synoptic on either screen displays the status of the plant components and also the sensor data required by ISO 7396-1 standard.




Main advantages :

- Fully automatic on-demand production
- Low operating costs
- Modular design
- Data logging and real time trend of events
- Remote monitoring capability

To assist operating and service, the generators are delivered with their full documentation. Multi-level secured access permits a supervisory control of several plants at once. A visual reminder of the maintenance times is also included on the Procom features.

For the commissioning of your plant, Mil's offers comprehensive start-up contracts where a technical specialist performs final arrangement and trains the end users.

VPSA oxygen generators



PRO₂XY[®] VAC - 5 bar

VPSA technology


- Totally oil-less process
- Oxygen purity 95 ± 1% or 93 ± 1%
- VPSA technology with low energy consumption
- Quick starting process, reduced time of 60% comparing to PSA
- Automatic calibration of the oxygen sensor(s) (option)
- Reduced maintenance
- Touch PLC control PROCOM3 developed according to EN62304
- Complies with EN 7396-1: 2016
- Environment-friendly



Characteristics

| PRO ₂ XY [®] VAC 5 bar version | Flow rate m ³ .h ⁻¹ 50 Hz - 20°C | | | Total installed power kW | | | Average power consumption - kW | | |
|---|---|------|------|-----------------------------|------|------|-----------------------------------|------|------|
| | 90% | 93 % | 95 % | 90% | 93 % | 95 % | 90% | 93 % | 95 % |
| PRO ₂ XY [®] VAC 27 | 34,6 | 30 | 27 | 23,7 | 23,7 | 23,7 | 16,7 | 16,3 | 15,7 |
| PRO ₂ XY [®] VAC 40 | 57 | 50 | 42 | 38,8 | 38,8 | 34,8 | 27,4 | 26,3 | 24,9 |

Process sequence



Legend
 Pneumatic



- ① Oil free low pressure compressor
- ② VPSA oxygen generator at 0,5 bar
- ③ Oil free vacuum pump
- ④ Buffer tank at 0,5 bar of oxygen
- ⑤ 5 bar O₂ compressor
- ⑥ 5 bar storage tank of oxygen
- ⑦ 5 bar hospital network

- HIGH PRESSURE OPTION
- ⑧ 200 bar - O₂ high pressure booster
 - ⑨ Rack of HP cylinders or cylinder manifolds
 - ⑩ High pressure filling skid for mobile cylinders

Power supply

3-phase, 400V+N+E / 50Hz



Range complying with
93/42/CEE directive



Range complying with
2014/68/UE directive
Réf : 520883-01-EN 01/2025

VPSA oxygen generators

Dimensions (mm)

PRO₂XY® VAC 27



PRO₂XY® VAC 40



| PRO ₂ XY® VAC 5 bar | PRO ₂ XY® VAC 27 | | PRO ₂ XY® VAC 40 | |
|-----------------------------------|-----------------------------|-------------|-----------------------------|-------------|
| | Dimensions (mm) | Weight (kg) | Dimensions (mm) | Weight (kg) |
| ① SRP | 1130 x 780 x 1200 | 300 | 1130 x 780 x 1200 | 325 |
| ② PAV | 1250 x 970 x 1300 | 425 | 1380 x 970 x 1300 | 475 |
| ③ Generator | 1525 x 1735 x 2190 | 1500 | 1710 x 1835 x 2250 | 2000 |

Compression module PRO₂XY® VAC 27



Compression module PRO₂XY® VAC 40



| PRO ₂ XY® VAC 5 bar | PRO ₂ XY® VAC 27 - 93/95 % | | PRO ₂ XY® VAC 40 - 93 % | | PRO ₂ XY® VAC 40 - 95 % | |
|-----------------------------------|---------------------------------------|-------------|------------------------------------|-------------|------------------------------------|-------------|
| | Dimensions (mm) | Weight (kg) | Dimensions (mm) | Weight (kg) | Dimensions (mm) | Weight (kg) |
| ⑤ Compression module | 1276 x 540 x 1273 | 170 | 1279 x 1115 x 1475 | 300 | 1276 x 540 x 1273 | 220 |
| ④ RT tank | 1200 x 880 x 2380 1000 l | 220 | 1490 x 1210 x 2520 2000 l | 395 | 1490 x 1210 x 2520 2000 l | 395 |
| ⑥ RP O ₂ tank | 1000 x 880 x 2380 1000 l | 220 | 1400 x 1210 x 2520 2000 l | 395 | 1400 x 1210 x 2520 2000 l | 395 |

Plant references

| PRO ₂ XY® VAC - 5 bar | | PRO ₂ XY® VAC 27 | | PRO ₂ XY® VAC 40 | |
|----------------------------------|--------------|-----------------------------|--------|-----------------------------|--------|
| | | 93 % | 95 % | 93 % | 95 % |
| PRO ₂ XY® VAC | Single line | 924972 | 924674 | 924976 | 924682 |
| Final filtration | Simple | 824211 | 824211 | 824217 | 824217 |
| PRO ₂ XY® VAC | Multi-line A | 924973 | 924675 | 924977 | 924683 |
| PRO ₂ XY® VAC | Multi-line B | 924974 | 924676 | 924978 | 924684 |
| PRO ₂ XY® VAC | Multi-line C | 924975 | 924677 | 924979 | 924685 |
| Final filtration | Duplex | 824212 | 824212 | 824218 | 824218 |

Accessories and option references

| | | | | |
|--|-----------------------|--------|-----------------------------------|--------|
| O ₂ paramagnetic analyzer | manual calibration | 724565 | O ₂ ambient analyzer | 622709 |
| O ₂ paramagnetic analyzer | automatic calibration | 724566 | O ₂ hygrometry sensor | 823735 |
| Redundancy box for O ₂ analyzer (EN62304) | | 724564 | CO KIT analyzer | 724912 |
| Redundancy box for HP (EN62304) | | 724752 | CO / CO ₂ KIT analyzer | 724914 |
| 2 nd O ₂ paramagnetic analyzer | manual calibration | 724627 | O ₂ flowmeter G1/2" | 624683 |
| 2 nd O ₂ paramagnetic analyzer | automatic calibration | 724628 | 24VCC power supply | 625004 |
| O ₂ zirconium analyzer | zirconium | 625616 | Energy monitoring - 100A | 625090 |

PRO₂XY[®] VAC - 12 bar

VPSA technology


- Totally oil-less process
- Oxygen purity 95 ± 1% or 93 ± 1%
- VPSA technology with low energy consumption
- Quick starting process, reduced time of 60% comparing to PSA
- Automatic calibration of the oxygen sensor(s) (option)
- Reduced maintenance
- Touch PLC control PROCOM3 developed according to EN62304
- Complies with EN 7396-1: 2016
- Environment-friendly



Characteristics

| PRO ₂ XY [®] VAC 12 bar version | Flow rate *m ³ .h ⁻¹ 50 Hz - 20°C | | | Total installed power kW | | | Average power consumption - kW | | |
|--|--|------|------|-----------------------------|------|------|-----------------------------------|------|------|
| | 90% | 93 % | 95 % | 90% | 93 % | 95 % | 90% | 93 % | 95 % |
| PRO ₂ XY [®] VAC 27 | 34,6 | 30 | 27 | 26,7 | 26 | 26 | 17,7 | 17,1 | 16,6 |
| PRO ₂ XY [®] VAC 40 | 57 | 50 | 42 | 43,3 | 42,6 | 38,6 | 29,8 | 28,3 | 27,1 |

Process sequence



- ① Oil free low pressure compressor
- ② VPSA oxygen generator at 0,5 bar
- ③ Oil free vacuum pump
- ④ Buffer tank at 0,5 bar of oxygen
- ⑤ 5 bar O₂ compressor
- ⑥ 5 bar storage tank of oxygen
- ⑦ 5 bar hospital network
- ⑧ 12 bar O₂ booster
- ⑨ 12 bar storage tank of oxygen
- ⑩ 12 bar hospital network

- HIGH PRESSURE OPTION
- ⑪ 200 bar - O₂ high pressure booster
 - ⑫ Rack of HP cylinders or cylinder manifolds
 - ⑬ High pressure filling skid for mobile cylinders

Power supply

3-phase, 400V+N+E / 50Hz



Range complying with
93/42/CEE directive



Range complying with
2014/68/UE directive
Réf : 520883-01-EN 01/2025

VPSA oxygen generators

Dimensions (mm)

PRO₂XY® VAC 27



PRO₂XY® VAC 40



| PRO ₂ XY® VAC 12 bar | PRO ₂ XY® VAC 27 | | PRO ₂ XY® VAC 40 | |
|------------------------------------|-----------------------------|-------------|-----------------------------|-------------|
| | Dimensions (mm) | Weight (kg) | Dimensions (mm) | Weight (kg) |
| ① SRP | 1130 x 780 x 1200 | 300 | 1130 x 780 x 1200 | 325 |
| ② PAV | 1200 x 970 x 1300 | 425 | 1380 x 970 x 1300 | 475 |
| ③ Generator | 1525 x 1735 x 2190 | 1500 | 1710 x 1835 x 2250 | 2000 |

Compression module PRO₂XY® VAC 27



Compression module PRO₂XY® VAC 40



| PRO ₂ XY® VAC 12 bar | PRO ₂ XY® VAC 27 - 93/95 % | | PRO ₂ XY® VAC 40 - 93 % | | PRO ₂ XY® VAC 40 - 95 % | |
|------------------------------------|---------------------------------------|-------------|------------------------------------|-------------|------------------------------------|-------------|
| | Dimensions (mm) | Weight (kg) | Dimensions (mm) | Weight (kg) | Dimensions (mm) | Weight (kg) |
| ⑤ Compression module | 1279 x 1060 x 1273 | 280 | 1279 x 1115 x 1620 | 410 | 1279 x 1060 x 1620 | 330 |
| ④ RT tank | 1200 x 880 x 2380 1000 l | 220 | 1490 x 1210 x 2520 2000 l | 395 | 1490 x 1210 x 2520 2000 l | 395 |
| ⑥ RP O ₂ tank | 1000 x 880 x 2380 1000 l | 220 | 1400 x 1210 x 2520 2000 l | 395 | 1400 x 1210 x 2520 2000 l | 395 |
| ⑦ R 10 tank | 970 x 880 x 2380 1000 l | 220 | 1410 x 1210 x 2520 2000 l | 395 | 1410 x 1210 x 2520 2000 l | 395 |

Plant references

| PRO ₂ XY® VAC - 12 bar | | PRO ₂ XY® VAC 27 | | PRO ₂ XY® VAC 40 | |
|-----------------------------------|--------------|-----------------------------|--------|-----------------------------|--------|
| | | 93 % | 95 % | 93 % | 95 % |
| PRO ₂ XY® VAC | Single line | 924980 | 924678 | 924984 | 924686 |
| Pressure release line | Simple | 824214 | 824214 | 824214 | 824214 |
| PRO ₂ XY® VAC | Multi-line A | 924981 | 924679 | 924985 | 924687 |
| PRO ₂ XY® VAC | Multi-line B | 924982 | 924680 | 924986 | 924688 |
| PRO ₂ XY® VAC | Multi-line C | 924983 | 924681 | 924987 | 924689 |
| Pressure release line | Duplex | 824215 | 824215 | 824215 | 824215 |

Accessories and option references

| | | | | |
|--|-----------------------|--------|--|--------|
| O ₂ paramagnetic analyzer | manual calibration | 724565 | O ₂ ambient analyzer | 622709 |
| O ₂ paramagnetic analyzer | automatic calibration | 724566 | O ₂ hygrometry monitoring kit | 823735 |
| Redundancy box for O ₂ analyzer (EN62304) | | 724564 | CO KIT analyzer | 724912 |
| Redundancy box for HP (EN62304) | | 724752 | CO / CO ₂ KIT analyzer | 724914 |
| 2 nd O ₂ paramagnetic analyzer | manual calibration | 724627 | O ₂ flowmeter G1/2" | 624683 |
| 2 nd O ₂ paramagnetic analyzer | automatic calibration | 724628 | Power supply 24VCC | 625004 |
| O ₂ probe analyzer kit | zirconium | 625616 | Energy monitoring - 100A | 625090 |